

## SPECIFICATION DATA


1-on-1 Conversations • TV Listening • Listening in the Car • Small Groups and much more

# PKT D1 Pocketalker® Ultra Personal Amplifier



The Pocketalker<sup>®</sup> Ultra amplifies sound clearly and easily for better understanding. Up to 200 hours of battery life ensures long-lasting performance. Each Pocketalker Ultra includes a plugin sensitive microphone and TV listening extension cord for the TV listener. Use the Pocketalker Ultra with a variety of available earphones and headphones of your choice, or a neckloop and a telecoil-equipped hearing aid.


## **PKT D1 Personal Amplifier:**




### Acoustic Specifications: Typical values (ANSI S3.22-1987)

#### Lightweight headphone, HED 021: (with 6cc coupler)


Test	Tone High	Tone Mid	Tone Low
SSPL90	124.3 dB	124.8 dB	124.8 dB
HF Ave (SSPL90)	120.7 dB	120.6 dB	120.5 dB
HF Ave full on gain (50 dB)	46.9 dB	44.9 dB	41.0 dB
Ref test gain	43.7 dB	43.6 dB	41.5 dB






#### Wide range headphone, HED 026: (with 6cc coupler)


Test	Tone High	Tone Mid	Tone Low
SSPL90	121.4 dB	125.3 dB	124.8 dB
HF Ave (SSPL90)	119.5 dB	120.7 dB	120.9 dB
HF Ave full on gain (50 dB)	45.6 dB	44.8 dB	41.2 dB
Ref test gain	42.5 dB	43.7 dB	41.7 dB






#### Mini earphone single, EAR 013: (with 2cc coupler)

Test	Tone High	Tone Mid	Tone Low
SSPL90	123.8 dB	124.0 dB	123.2 dB
HF Ave (SSPL90)	118.7 dB	118.6 dB	118.6 dB
HF Ave full on gain (50 dB)	44.7 dB	44.2 dB	38.9 dB
Ref test gain	41.7 dB	41.6 dB	39.2 dB





Note: The Pocketalker® Ultra is not sold, advertised or intended for use as a personal hearing aid. Williams AV recommends that any individual experiencing hearing loss be examined by a physician/audiologist to determine the cause before using this device.

#### **Bid Specs**

The Williams AV Pocketalker Ultra (model PKT D1) shall be housed in an ABS plastic, impact-resistant, silver case with a removable battery door. The dimensions shall be 90mm X 55mm X 20mm.

The Pocketalker Ultra shall operate on 2x AAA batteries and operate for up to 200 hours continually with disposable alkaline batteries.

The Pocketalker Ultra shall have a maximum 130dB output with EAR 013 earphone and have a non-adjustable AGC (Limiter) and external rotary adjust tone control. The Pocketalker Ultra shall have a combination ON/OFF rotary volume control and a red power LED that will flash to signal when the battery level is low. The Pocketalker Ultra shall have 3.5mm mono microphone and earphone jacks.

The Pocketalker Ultra shall have a removable belt clip for handsfree operation and a lanyard connection at the bottom of the case. The Pocketalker Ultra shall have CE approval and be compliant with RoHS and WEEE regulations.

The Pocketalker Ultra shall have a five-year warranty excluding earphones, headphones, batteries, chargers, belt clip and other accessories.

The Pocketalker Ultra shall be the Williams AV model PKT D1.

#### **Domestic Sales**

Williams AV 10300 Valley View Rd Eden Prairie, MN 55344 Phone: 800-328-6190 / 952-943-2252 Fax: 952-943-2174 Email: info@williamsav.com Web: www.williamsav.com

#### **International Sales**

International Sales Department Williams AV 10300 Valley View Rd Eden Prairie, MN 55344 USA Phone: +1 952 943 2252 Fax: +1 952 943 2174 Email: info@williamsav.com Web: www.williamsav.com

